


Dialysate compartment

Outer belt, contains channel housing for CAPD dialysis components and tubing


Outer belt channel section for Stay-Safe® organizer

**Figure 1:**  
User wearing the MobileExchange Pack, front view




**Figure 2: User wearing the MobileExchange Pack, rear view**


Tubing channel and organizer housing. All dialysis components and tubing are accessible to the user only



**Figure 3:**  
User wearing the MobileExchange Pack,  
front view, belt open




**Figure 4:**  
**MobileExchange Pack,**  
**front view, features breakdown**





**Figure 5:**  
**MobileExchange Pack, front view, open,**  
**with stay-safe® CAPD system placement**

When empty and not in use, the drain bag is rolled up and tucked inside the belt "channel"



DRAIN BAG




Velcroed straps and loops keep Stay-Safe organizer and tubing securely in place

**Figure 6:**  
**MobileExchange Pack, front view, open,**  
**drain bag placement**




Pouch for drain bag, available for use during drain

**Figure 7a:**  
**MobileExchange Pack,**  
**drain bag pouch accessory, stored**



**Figure 7b:**  
**MobileExchange Pack,**  
**drain bag pouch accessory, in use**

Drain bag pouch hangs from loops at bottom of outer belt, and provides mobility and privacy during drain.



**Figure 8:**  
**MobileExchange Pack, fill/drain**  
**fluid routing map**